



Conforms to: IS / ISO 9170 - 1 : 2017, IS / ISO 9170 - 2 : 2008

DIN 13620 - 2, BS 5682 : 2015 & CGA V-5-2019

Certifications: ISO 13485:2016 – Medical devices – Quality Management Systems ISO 9001:2015 - Quality Management Systems

Special Features

- Stable Connection
- Quick to connect and disconnect
- 3 delta finish for long O-ring life of the outlet
- Indexed to specific gas service
- · Laser engraved identification of gas service
- · Customizable patient end connection
- Right angled (90 degree) patient end connection on request
- MRI compatible
- Made from AISI 303 Steel

Probes are crucial for gas specific and stable connection. These are used to connect medical gases and vacuum devices to outlets / terminal units of MGPS distribution network.

Manufactured in wide range of Quick connect and DISS types compatible to various Outlet / Terminal Unit standards, maintaining gas specific keying safety.

Patient end connection to medical devices are available in 6 mm OD Hose barb, 1/8" BSPT, 3/8" BSP etc. However, customized patient end connection, angled connection and various end fittings are also available on request.

Probes are available for O_2 , N_2O , O_2 / N_2O (Entonox), CO_2 , Vacuum, Medical Air -AR4, Instrument Air -AR7, Nitrogen, Helium/oxygen mixtures (Heliox), AGS, Triple mix gases (Tri-Gas: CO_2 , $O_2 & N_2$) and gas mixtures classified as medical gases.

Ordering Information

British Standard	Quick Conne	ect Type	Keying Style - Gas Specific Diameter Indexing								
(BS 5682:2015)		End Gas Services									
	Probe	Connections	O ₂	N₂O	AR 4	AR 7	VAC	AGSS*	CO ₂ *		
		6mm OD Hose barb	6013101	6013401	6013301	6013501	6013201	6013701	6013601		
		3/8" BSP	6013103	6013403	6013303	6013503	6013203	xx	6013603		

^{*}CO₂ & AGSS Outlets include a NIST style connection as there is no defined probe size in British Standard.

	Quick Connect T	уре	Keying Style - Gas Specific				ter Indexin	g			
	Probe	End	Gas Services								
DIN Standard	Probe	Connections	O ₂	N₂O	AR 4	AR 7	VAC	CO ₂			
13260-2		6mm OD Hose barb	6018101	6018401	6018301	6018501	6018201	6018601			
•		3/8" BSP	6018103	6018403	6018303	6018503	6018203	6018603			
7.		1/8" BSPT	6018102	6018402	6018302	6018502	6018202	6018602			

DISS	Screw Connect	t Type	Keying Style - Gas Specific Diameter Index Safety System								
CGA V-5-2019		End	End Gas Services								
	Probe	Connections	O ₂	N₂O	AR 4	AR 7	VAC	AGSS	CO ₂		
Miller Sycology		6mm OD Hose barb	6016101	6016401	6016301	6016501	6016201	6016701	6016601		
		1/8" BSPT	6016102	6016402	6016302	6016502	6016202	XX	6016602		

_ '

Ohmeda (Diamond)	Quick Connect T	уре	Keying Style - Gas Specific Notch Index					
Campatible	Probe	End			Gas S	ervices		
		Connections	O ₂	N₂O	AR 4	AR 7	VAC	CO ₂
		6mm OD Hose barb	6014101	6014401	6014301	On Request	6014201	On Request
		1/8" BSPT	6014102	6014402	6014302	On Request	6014202	On Request

	Quick Connec	t Type	Keying Style - Gas Specific Geometric Indexir					ic Indexin	g			
Puritan Bennett (PB) Compatible		End	Gas Services									
	Probe	Connections	O ₂	O ₂ N ₂ O AR 4 AR 7 V	VAC	AGSS	CO ₂					
NITROUS ONICE		6mm OD Hose barb	6011101	6011401	6011301	6011501	6011201	6011701	6011601			
		1/8" BSPT	6011102	6011402	6011302	6011502	6011202	6011702	6011602			
		3/8" BSP	6011103	On Request	On Request	On Request	6011203	On Request	On Request			

Parkodex (Emerald)	Quick Connect T	Туре	Keying Style - Gas Specific Pin Index Safety System (P					
Compatible		End			Gas Se	ervices		
	Probe	Connections	O ₂	N₂O	AR 4	AR 7	VAC	CO ₂
		6mm OD Hose barb	6012101	6012401	6012301	6012501	6012201	6012601
		3/8" BSP	6012103	6012403	6012303	6012503	6012203	6012603

Anaesthetic Gas Scavenging Probes





Type - BS Outlet

Connection: 13/16 x 13 TPI Screw Connect Probe Code: 6013701





Type - DISS Outlet

CGA V5 Connection Number 2220

Screw Connect Probe Code: 6016701

Patient End Conversion Connectors





Code: 6019001 Code: 6019002

DISCLAIMER: Aktiv Probes / Connectors / Adapters meet and exceed the requirements of relevant standards at the time of manufacture. However piping source, capacity, sizing and restrictions may prevent from attaining the stipulated value.



AKTIV TECHNOLOGIES

Corporate Office : Plot No. 691, Sector 69, IMT, Faridabad, Haryana 121004, INDIA Tel: +91-7291838253

Email: aktivtechnologies@gmail.com, sanjay.aktiv@gmail.com

Website: www.aktivtechnologies.com



_IT / 601A Rev 02 12/2