



## Features


- Reusable, efficient and hygienic
  - Mechanical overfill cut-off protection
  - Made of shock-resistant polycarbonate
  - Suitable for normal disinfection & sterilization
  - Autoclavable at 121°C for 15 minutes
  - Leak proof seal to prevent vacuum loss
  - Clearly distinctive patient and vacuum ports
  - Clear, easy to read, highly visible graduated markings
  - Easy to set-up with locking clasp to secure jar on to the lid
- Provide safe liquid medical waste disposal
  - Suitable for the collection of small and large volumes of liquids
  - Enables hygienic collection of aspirated secretion, suctioned fluids and other organic fluids safely

**Conforms to:** IS / ISO 10079 - 3 : 2014, IS / ISO 7396 - 1 : 2016 + A1 : 2017  
HTM 02-01 Part A - 2006

**Certifications:** ISO 9001:2015 - Quality Management Systems  
ISO 13485:2016 – Medical devices – Quality Management Systems

# Technical Specifications

Conforms to: IS / ISO 10079 - 3 : 2014, IS / ISO 7396 - 1 : 2016 + A1 : 2017, HTM 02-01 Part A - 2006



Vacuity Collection Jar Variants	Capacity (Litre)		Dimensions				Ordering Code
	Jar	Fill	A	B	C	D	
J800	1	0.8	127	105	168	195	6022140
J800 DISS	1	0.8	127	105	168	217	6022141
J1000	1.75	1.0	160	140	150	181	6022150
J1000 DISS	1.75	1.0	160	140	150	203	6022151
J1750	2.25	1.75	160	140	214	244	6022160
J1750 DISS	2.25	1.75	160	140	214	266	6022161
J3500	4	3.5	160	140	363	394	6022170
J3500 DISS	4	3.5	160	140	363	416	6022171



Clean, Disinfect / Sterilize the jar before and after each use.



**AKTIV TECHNOLOGIES**

**Corporate Office :** Plot No. 691, Sector 69, IMT, Faridabad, Haryana 121004, INDIA  
Tel: +91-7291838253  
Email: [aktivtechnologies@gmail.com](mailto:aktivtechnologies@gmail.com), [sanjay.aktiv@gmail.com](mailto:sanjay.aktiv@gmail.com)  
Website : [www.aktivtechnologies.com](http://www.aktivtechnologies.com)



Consult Instructions for Use

LT / 602C Rev 00 12/21

**Disclaimer:** Due to our policy of continuous product improvement, Aktiv Technologies reserve the right to alter dimensions and/or specification without notice although every effort will be made to advise of any such modifications.